



Email : sales@laboratoryinstrumentindia.com

Product Name :
Conservation Of Linear Momentum

Product Code :
LBNY-0005-13300010



Description :

Versatile unit allowing the study of Conservation of Linear momentum
To comprise of two trolleys which runs on a precision air track which can be inclined to the horizontal
Trolleys projected along air track using catapult mechanism
Trolley mass to be adjusted by adding weights
To be a versatile apparatus which covers other kinetic principles such as motion, acceleration, Newton's laws of motion, velocity
Comprehensive technical manual
Set of slotted weights included
To be supplied with timing gates and timer.

Technical Specification :

Conservation Of Linear Momentum



Laboratory Instrument India