



Email : sales@laboratoryinstrumentindia.com

Product Name :
Blood Pressure Sensor

Product Code :
LBNY-0004-780044



Description :

Blood Pressure Sensor

Technical Specification :

Blood Pressure sensor:

Measures heart rate and systolic and diastolic arterial blood pressure (mmHG)

Sensor used:

To determine the effects of exercise on blood pressure and heart rate

To explore the effects of body position on blood pressure and heart rate

Include:

Blood pressure sensor

Standard arm cuff

Bladder and pressure release valve

Supplied with extension cable.



Laboratory Instrument India

