



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Blood Pressure Sensor

**Product Code :**  
LBNY-0004-780044



**Description :**

Blood Pressure Sensor

**Technical Specification :**

Blood Pressure sensor:

Measures heart rate and systolic and diastolic arterial blood pressure (mmHG)

Sensor used:

To determine the effects of exercise on blood pressure and heart rate

To explore the effects of body position on blood pressure and heart rate

Include:

Blood pressure sensor

Standard arm cuff

Bladder and pressure release valve

Supplied with extension cable.



**Laboratory Instrument India**

