




Email : sales@laboratoryinstrumentindia.com

Product Name :
Shearing Force Apparatus

Product Code :
LBNY-0005-13300018



Description :

The bench mounted apparatus consists of a 'shear beam', 2 supports and 3 load hangers. The 'shear beam' of two parts, joined at a low friction joint. Internal forces created, and bending moment resisted by an underslung adjustable tension spring: an overhung spring balance provides the vertical shearing force. Demonstration of the action of shear in a beam; measurement of shear force at a section of a loaded beam; verify equilibrium of vertical forces. The load hangers mountable in virtually any position along the beam. An Instruction manual for student and lecturer provided. Set of weights.

Technical Specification :

Shearing Force Apparatus



Laboratory Instrument India