




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Bell Crank Lever

**Product Code :**  
LBNY-0005-13300021



**Description :**

The lever arm is with scale and grooved at every 5 cm. complete with one sliding weight of 1 kg. A wooden right angled lever with fulcrum at the angle, the horizontal arm is of 75 cm. and vertical arm is of 20 cm., a spring balance with adjustable wing nut is attached to the shorter arm.

**Technical Specification :**

Bell Crank Lever



**Laboratory Instrument India**