



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bending Moment Apparatus

Product Code :
LBNY-0005-13300022



Description :

The experimental beam is in two parts, joined by a pair of low friction ball bearings. An underslung spring balance provides a resisting moment and measures the section bending moment.
Demonstration of the action of moment of resistance in a beam; measurement of bending moment at a section of a loaded beam; verify equilibrium of forces and moments.
The load hangers mountable in virtually any position along the beam.
An Instruction manual for student and lecturer provided
To be a bench mounted apparatus consisting of an 'experimental beam', 3 weight hangers, 2 supports and a spring balance
Set of calibrated weights and tools supplied.

Technical Specification :

Bending Moment Apparatus



Laboratory Instrument India