



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bending Moment Apparatus

**Product Code :**  
LBNY-0005-13300022



**Description :**

The experimental beam is in two parts, joined by a pair of low friction ball bearings. An underslung spring balance provides a resisting moment and measures the section bending moment.  
Demonstration of the action of moment of resistance in a beam; measurement of bending moment at a section of a loaded beam; verify equilibrium of forces and moments.  
The load hangers mountable in virtually any position along the beam.  
An Instruction manual for student and lecturer provided  
To be a bench mounted apparatus consisting of an 'experimental beam', 3 weight hangers, 2 supports and a spring balance  
Set of calibrated weights and tools supplied.

**Technical Specification :**

Bending Moment Apparatus



**Laboratory Instrument India**