



Email : sales@laboratoryinstrumentindia.com

Product Name :
Toxicity Apparatus

Product Code :
LBNY-0005-1350003



Description :

Features:

An internally mounted stirring fan for rapid mixing of combustion products
12 sampling ports on the side of the chamber
The gas burner has a spark ignition system which automatically reignites should the flame extinguish
A forced-air extraction system for evacuating the chamber after a test
A separate control unit houses the flowmeters, timer, methane and air controls
The combustion chamber has welded seams and a volume of 0.7m³.

Technical Specification :

Toxicity Apparatus



Laboratory Instrument India