




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Light Sensor

**Product Code :**  
LBNY-0004-780045



**Description :**

Light Sensor

**Technical Specification :**

Light Sensor:

The High Sensitivity Light Sensor is designed to be used with a compatible interface to make measurements of relative light intensity. It is designed to perform visible light studies from low-intensity spectral studies to daylight. Built-in automatic variable oversampling reduces noise.

Measure in bright sunlight to dim candlelight, and everything in between. Light Sensor has three ranges to accommodate a wide variety of student light measurements. Switch ranges with the touch of a button and without losing calibration. No other Light Sensor combines such a wide range with such high-sensitivity.

Photodiode spectral response from 320 to 1100nm

Three switches selectable gain levels: 1X, 10X and 100X

Supplied with extension cable

Provide necessary accessories.



**Laboratory Instrument India**

