



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Set Of Sectioned Components For Pneumatic Controls

**Product Code :**  
LBNY-0005-1370009



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The set includes some of the most important pneumatic components, used in industrial control applications. Special explanatory technical plates show each sectioned component explaining its operation with notes, drawings and technical data. The use of contrasting colors according to a pre-established code simplifies the informative process as well as the learning one. In sectionizing each component the best care has been taken to evidence all the interesting constructive detailed features.

**Technical Specification :**

Set Of Sectioned Components For Pneumatic Controls



**Laboratory Instrument India**