



Email : sales@laboratoryinstrumentindia.com

Product Name :
Light Sensor High Sensitivity

Product Code :
LBNY-0004-780046



Description :

Light Sensor High Sensitivity

Technical Specification :

Light Sensor (high sensitivity)
Photodiode spectral response from 320 to 1100nm
Three switches selectable gain levels: 1X, 10X and 100X
Light intensity levels: 5Lux @ 1X, 0.5Lux @ 10X and 0.05Lux@ 100X gain
Provide necessary accessories.



Laboratory Instrument India