



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Light Sensor High Sensitivity

**Product Code :**  
LBNY-0004-780046



**Description :**

Light Sensor High Sensitivity

**Technical Specification :**

Light Sensor (high sensitivity)  
Photodiode spectral response from 320 to 1100nm  
Three switches selectable gain levels: 1X, 10X and 100X  
Light intensity levels: 5Lux @ 1X, 0.5Lux @ 10X and 0.05Lux@ 100X gain  
Provide necessary accessories.



**Laboratory Instrument India**