



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Fluid Threads Visualizer

**Product Code :**  
LBNY-0005-13700011



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It consists of a transparent illuminated test chamber where the aerodynamic models are mounted and of a smoke generator. Fluid Threads Visualizer has been designed in order to visualize the fluid threads in a low-turbulence wind tunnel.

The Fluid Threads Visualizer unit is available in two versions:

Fluid threads visualizer with manual control

Fluid threads visualizer with motor control

**Technical Specification :**

Weight and dimensions of Fluid Threads Visualizer

Dimensions: 1400 x 550 x 500 h mm

Weight: 40 kg.

**Required services:**

Electric supply: 220V single-phase, 50/60 Hz; 0.5 kW.



**Laboratory Instrument India**

