



Email : sales@laboratoryinstrumentindia.com

Product Name :
Variable Speed Two Stage Compressor Study Unit With
Digital Torque Indicator

Product Code :
LBNY-0005-13700014



www.laboratoryinstrumentindia.com

Description :

The bench is trolley-mounted and has a front panel on which the control devices and all the instruments needed for the tests are located. The system includes two water coolers to cool the air generated by the compressor, complete with thermo-resistances for measuring the temperature of the air and water at the input and output. The Variable Speed Two Stage Compressor Study Unit With Digital Torque Indicator makes it possible to study the major characteristics of a twin-cylinder air compressor with two-stage compression.

Technical Specification :

Weights and dimensions:

Dimensions: 2100 x 950 x 2000h mm

Weight: 450 kg.

Required services:

Water supply: 600l/h, 2 - 3bar

Electrical power supply: 220 V single-phase, 50/60 Hz, 3 kW.



Laboratory Instrument India

