



Email : sales@laboratoryinstrumentindia.com

Product Name :
Compressors Study Unit Reciprocating And Rotative

Product Code :
LBNY-0005-13700016



Description :

It is supplied complete and self-sufficient, needing only the connection to the electrical feeding to the water network. The Compressors Study Unit Reciprocating And Rotative allows the study of the working characteristics of reciprocating and rotative compressors.

Technical Specification :

Required services:

Electric supply: 220/240 V, 50/60 Hz, single-phase, 4 KW

Weight and dimension:

Module for compressed air tests: 1500 x 700 x 1800 h mm - 450 kg

Compressors module: 1000 x 700 x 2300 h mm - 300 kg.



Laboratory Instrument India