




Email : sales@laboratoryinstrumentindia.com

Product Name :
Double Stage Compressor Study Unit

Product Code :
LBNY-0005-13700018



Description :

The system includes two water refrigerators to cool the air generated by the compressor, complete with thermo-resistances for measuring the temperature of the air and water at the input and output. The Double Stage Compressor Study Unit makes it possible to study the main characteristics of a twin-cylinder air compressor with two-stage compression. The bench is trolley-mounted and has a front panel on which the control devices and all the instruments needed for the tests are located.

Technical Specification :

Required services:

Water supply: 600 l/h, 2 - 3bar

Electrical power supply: 380 V three-phase + neutral, 50/60 Hz, 3 kW

Weights and dimensions:

Dimensions: 2100 x 950 x 2000h mm

Weight: 415 kg.



Laboratory Instrument India

