




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Double Stage Compressor Study Unit

**Product Code :**  
LBNY-0005-13700018



**Description :**

The system includes two water refrigerators to cool the air generated by the compressor, complete with thermo-resistances for measuring the temperature of the air and water at the input and output. The Double Stage Compressor Study Unit makes it possible to study the main characteristics of a twin-cylinder air compressor with two-stage compression. The bench is trolley-mounted and has a front panel on which the control devices and all the instruments needed for the tests are located.

**Technical Specification :**

**Required services:**

Water supply: 600 l/h, 2 - 3bar

Electrical power supply: 380 V three-phase + neutral, 50/60 Hz, 3 kW

**Weights and dimensions:**

Dimensions: 2100 x 950 x 2000h mm

Weight: 415 kg.



**Laboratory Instrument India**

