



Email : sales@laboratoryinstrumentindia.com

Product Name :
Axial Fan Module

Product Code :
LBNY-0005-13700019



Description :

Pressure measurement: Electronic transducers and digital display
Power, speed and torque measurement:
From Universal Dynamometer
Nett dimensions and weight (assembled):
2300 mm x 1060 mm x 1530 mm and 200 kg
Packed volume and weight: 2.4 m³ and approximately 300 kg
Nominal Fan Diameter: 400 mm
Flow rate: 1.7 m³.s⁻¹

Technical Specification :

Axial Fan Module



Laboratory Instrument India