



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Axial Fan Module

**Product Code :**  
LBNY-0005-13700019



**Description :**

Pressure measurement: Electronic transducers and digital display  
Power, speed and torque measurement:  
From Universal Dynamometer  
Nett dimensions and weight (assembled):  
2300 mm x 1060 mm x 1530 mm and 200 kg  
Packed volume and weight: 2.4 m<sup>3</sup> and approximately 300 kg  
Nominal Fan Diameter: 400 mm  
Flow rate: 1.7 m<sup>3</sup>.s<sup>-1</sup>

**Technical Specification :**

Axial Fan Module



**Laboratory Instrument India**