



Email : sales@laboratoryinstrumentindia.com

Product Name :
Optical Dissolved Oxygen Sensor

Product Code :
LBNY-0004-780047



Description :

Optical Dissolved Oxygen Sensor

Technical Specification :

Optical Dissolved Oxygen Sensor:
Sealed body-electrode to allow direct submersion in aquatic environments.
Oxygen concentration range: 0 mg/L - 20mg/L min or 0-300% saturation
Response time: 90% in 25 sec
No stirring required
No warm up time
No calibration required
Built in pressure and temperture compensation
3m cable length
Can report of ambient pressure
Supplied with metal guard.



Laboratory Instrument India

