



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Optical Dissolved Oxygen Sensor

**Product Code :**  
LBNY-0004-780047



**Description :**

Optical Dissolved Oxygen Sensor

**Technical Specification :**

Optical Dissolved Oxygen Sensor:  
Sealed body-electrode to allow direct submersion in aquatic environments.  
Oxygen concentration range: 0 mg/L - 20mg/L min or 0-300% saturation  
Response time: 90% in 25 sec  
No stirring required  
No warm up time  
No calibration required  
Built in pressure and temperture compensation  
3m cable length  
Can report of ambient pressure  
Supplied with metal guard.



**Laboratory Instrument India**

