



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifugal Compressor Module

Product Code :
LBNY-0005-13700020



Description :

Experiments

Determination of the specific speed of the fan Essential Base Unit
Performance of a centrifugal fan
Universal Dynamometer
Variation of fan performance with type of impeller
Variation of fan performance with speed
Non-dimensional performance curves.

Recommended Ancillaries:

Pipe Flow and Nozzle Kit
VDAS-F (frame mounted version of the Versatile Data Acquisition System).

Technical Specification :

Essential Services

Electrical supply (for the Universal Dynamometer):
Phase to phase supply - 220 VAC, 60 Hz at 20 A
Phase to neutral supply - 230 VAC, 50 Hz at 20 A.

Laboratory Instrument India

