



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Oxygen Gas Digital Sensor

**Product Code :**  
LBNY-0004-780048



**Description :**

Oxygen Gas Digital Sensor

**Technical Specification :**

Oxygen gas digital Sensor:  
Measures gas concentration in open or closed system  
Range between 0 to 100%  
Accuracy: +/- 1%  
Resolution: 0.025%  
Supplied with extension cable.



**Laboratory Instrument India**