



Email : sales@laboratoryinstrumentindia.com

Product Name :
Surface Tension Apparatus

Product Code :
LBNY-0005-1380009



Description :

A horizontal ring is suspended on load cell and immersed in a liquid in a dish which is put on an adjustable height support stand. Surface Tension Apparatus is designed for measuring the surface tension of liquids. The force can be displayed on a force meter.

Technical Specification :

Frame : Aluminum profile with rubber footings.
Load cell : Two directions with hook.
Measuring ring : Plastic rings with suspension chords.
Support stand : Adjustable height.
Glass dish : 1 ea.
Power supply : 4 x 1.5 V AA battery or Adaptor.
Force measurement : Digital display in Newton or gram.
a. Measuring range : 200 mN/20g.
b. Zero compensation : Manual.



Laboratory Instrument India

