



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Surface Tension Apparatus

**Product Code :**  
LBNY-0005-1380009



**Description :**

A horizontal ring is suspended on load cell and immersed in a liquid in a dish which is put on an adjustable height support stand. Surface Tension Apparatus is designed for measuring the surface tension of liquids. The force can be displayed on a force meter.

**Technical Specification :**

Frame : Aluminum profile with rubber footings.  
Load cell : Two directions with hook.  
Measuring ring : Plastic rings with suspension chords.  
Support stand : Adjustable height.  
Glass dish : 1 ea.  
Power supply : 4 x 1.5 V AA battery or Adaptor.  
Force measurement : Digital display in Newton or gram.  
a. Measuring range : 200 mN/20g.  
b. Zero compensation : Manual.



**Laboratory Instrument India**

