




Email : sales@laboratoryinstrumentindia.com

Product Name :
Two And Three Pinned Arch

Product Code :
LBNY-0005-13800014



Description :

Experiments:

Characteristics of a three-pinned arch

Appreciation of footing stability and economy

Operating Conditions

Relationship between applied loads and horizontal thrust produced from a simple determinate arched structure.

Operating environment: Laboratory

Operating relative humidity range: 80% at temperatures Storage temperature range: 25°C to +55°C (when packed for transport)

Operating temperature range: +5°C to +40°C.

Standard Features:

Made in accordance with the latest European Union directives.

Two And Three Pinned Arch

Supplied with lecturer guide and student guide.

Technical Specification :

Two And Three Pinned Arch

Laboratory Instrument India

