




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Two And Three Pinned Arch

**Product Code :**  
LBNY-0005-13800014



**Description :**

**Experiments:**

Characteristics of a three-pinned arch

Appreciation of footing stability and economy

Operating Conditions

Relationship between applied loads and horizontal thrust produced from a simple determinate arched structure.

Operating environment: Laboratory

Operating relative humidity range: 80% at temperatures Storage temperature range: 25°C to +55°C (when packed for transport)

Operating temperature range: +5°C to +40°C.

**Standard Features:**

Made in accordance with the latest European Union directives.

Two And Three Pinned Arch

Supplied with lecturer guide and student guide.

**Technical Specification :**

Two And Three Pinned Arch

---

# Laboratory Instrument India

