



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wireless Carbon Dioxide Gas Sensor

Product Code :
LBNY-0004-780049



Description :

Wireless Carbon Dioxide Gas Sensor

Technical Specification :

Wireless Carbon Dioxide Gas Sensor
Measures gas concentration in open or closed system
Resolution: 2 ppm
Range between 0 to 100000 ppm
Accuracy +/- 15% full scale
Can be also connected via USB
Supplied with waterproof sleeve (at least quantity 20).



Laboratory Instrument India