




Email : sales@laboratoryinstrumentindia.com

Product Name :
Bending Stress In A Beam

Product Code :
LBNY-0005-13800018



Description :

Study of Bending Stress In A Beam
Converting strains to stresses
Strain gauges
Second moment of area
The bending equation
The neutral axis.

Technical Specification :

Operating Conditions:

Operating environment: Laboratory
Operating relative humidity range: 80% at temperatures Storage temperature range: 25°C to +55°C (when packed for transport)
Operating temperature range: +5°C to +40°C.



Laboratory Instrument India

