




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bending Stress In A Beam

**Product Code :**  
LBNY-0005-13800018



**Description :**

Study of Bending Stress In A Beam  
Converting strains to stresses  
Strain gauges  
Second moment of area  
The bending equation  
The neutral axis.

**Technical Specification :**

**Operating Conditions:**

Operating environment: Laboratory  
Operating relative humidity range: 80% at temperatures Storage temperature range: 25°C to +55°C (when packed for transport)  
Operating temperature range: +5°C to +40°C.



**Laboratory Instrument India**

