




Email : sales@laboratoryinstrumentindia.com

Product Name :
Plastic Bending Of Beams Apparatus

Product Code :
LBNY-0005-13800019



Description :

Experiments:

Introduction to limit state design

Relationship between maximum loading and plastic hinge formation for a simply supported beam, a propped cantilever and a fixed beam

Relationship between load and deflection for beams loaded to the plastic condition

Introduction to form factor.

Technical Specification :

Technical Specification:-

Operating temperature range: +5°C to +40°C

Operating relative humidity range: 80% at temperatures Operating environment: Laboratory

Storage temperature range: 25°C to +55°C (when packed for transport).



Laboratory Instrument India