



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Unsymmetrical Bending And Shear Centre

**Product Code :**  
LBNY-0005-13800020



**Description :**

**Experiments:**

The relationship between the vertical and horizontal deflections and the principal moments of area of each section.

The shear centre of various asymmetrical sections.

Horizontal and vertical deflection of different asymmetrical sections at various angles.

Horizontal and vertical deflection of different asymmetrical sections under various loads.

**Technical Specification :**

**Operating Conditions:**

Operating environment: Laboratory

Operating relative humidity range: 80% at temperatures Storage temperature range: 25°C to +55°C (when packed for transport)

Operating temperature range: +5°C to +40°C.



**Laboratory Instrument India**

