




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Continuous And Indeterminate Beam

**Product Code :**  
LBNY-0005-13800021



**Description :**

**Experiments:**

Reactions and fixing moments of a fixed beam and a propped cantilever  
Reaction and fixing moment of a propped cantilever with a sinking support  
Reactions of a simply supported beam  
Reactions of a two-span continuous beam  
This equipment allows many possible experiment  
Relationship between load and deflection for beams and cantilevers  
Configurations, using a stiff (rigid) beam, or a significantly more flexible beam.

**Technical Specification :**

**Operating Conditions:**

Operating environment: Laboratory  
Operating temperature range: +5C to +40C  
Operating relative humidity range: 80% at temperatures Storage temperature range: 25C to +55C (when packed for transport).



**Laboratory Instrument India**

