




Email : sales@laboratoryinstrumentindia.com

Product Name :
Continuous And Indeterminate Beam

Product Code :
LBNY-0005-13800021



Description :

Experiments:

Reactions and fixing moments of a fixed beam and a propped cantilever
Reaction and fixing moment of a propped cantilever with a sinking support
Reactions of a simply supported beam
Reactions of a two-span continuous beam
This equipment allows many possible experiment
Relationship between load and deflection for beams and cantilevers
Configurations, using a stiff (rigid) beam, or a significantly more flexible beam.

Technical Specification :

Operating Conditions:

Operating environment: Laboratory
Operating temperature range: +5C to +40C
Operating relative humidity range: 80% at temperatures Storage temperature range: 25C to +55C (when packed for transport).



Laboratory Instrument India

