




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Automatic Data Acquisition Unit

**Product Code :**  
LBNY-0005-1390003



**Description :**

It allows data logging, analysis and extra virtual simulated experiments. It accepts inputs from a digital force display, a digital strain display, an angular sensor and digital deflection indicators. It converts these inputs into the correct data signals for the computer. The Automatic Data Acquisition Unit is an interface box with software that connects experiments from the Structures range to a suitable computer. The software can also simulate experiments which students can perform using the hardware, so they can compare simulated and real results. The software can then analyse the data and create tables and charts.

**Technical Specification :**

Automatic Data Acquisition Unit



**Laboratory Instrument India**