



Email : sales@laboratoryinstrumentindia.com

Product Name :
Analog Lab Unit

Product Code :
LBNY-0002-610008



Description :

Analog Lab Unit

Technical Specification :

High level, high quality digital-analog trainer
Replaceable 4 pin connector, easy to maintenance
Combines all essential function of analog and digital experiment
With removable breadboard, DC power supply, function generator, two pulse switches, 2 1/2 inch 8 ohm 0.5W speaker and etc.
Fixed DC: +5V, +12V, -12V
Maximum Current Output: 0.3A
Dual DC Power Supply: 0 ~ +25V and 0 ~ -25V
Adjustable AC Power Supply
Voltage Range: 2V-0V-2V, 12V-0V-12V, 15V-0V-15V
Maximum Current Output: 500mA
Function Generator:
Output Waveform: Sine, Square, Triangle and TTL
Output Frequency: up to 100KHz in five steps

Analog Lab Unit, Analog Lab Unit Equipments, Analog Lab Unit Tools, Analog Lab Unit Tool Kits, Analog Lab Unit Manufacturers, Analog Lab Unit Suppliers from India, China, Kenya



Laboratory Instrument India