



Email : sales@laboratoryinstrumentindia.com

Product Name :
Universal Frame and Stand

Product Code :
LBNY-0005-13900044



Description :

The floor mounting frame can be easily converted to a bench-mounted version when using the bench mounting feet arrangement which is supplied as standard. The internal working dimensions for mounting experiments are 1.22m x 0.76m. Manufactured from high quality, ergonomic aluminium profile, the Universal Frame and Stand comprises an upper dual frame, which creates a continuous mounting slot around all four sides of the frame thus allowing the individual experiment components to be interchanged quickly. Should specimen beams longer than 1.22metres be used, the working gap allows the beams to run through the frame and protrude outside the working area.

Technical Specification :

To be easily converted to bench mounted frame using parts supplied as standard
To allow attachment of Interface unit via profiled grooves of frame
Enables student project work
To be floor standing, twinned aluminium profile section frame
To have 5000N (5kN) capacity
To be supplied with full comprehensive instruction manual and tools
To mount all structures experiments quickly and easily
To have 1.22 x 0.76m clear internal work space.

Laboratory Instrument India

