



Email : sales@laboratoryinstrumentindia.com

Product Name :
Suspension Bridge Apparatus

Product Code :
LBNY-0005-13900051



Description :

The apparatus is equally suitable for student experiments or demonstrations. Model consists of a three-hinged deck suspended by ties from a light cable. This bench-mounted apparatus is used to investigate the structure under a wide range of load conditions. Cable can be readily adjusted to give different values of dip over the 1m span.

Technical Specification :

Suspension Bridge Apparatus



Laboratory Instrument India