



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Drop Counter Sensor

**Product Code :**  
LBNY-0004-780055



**Description :**

Drop Counter Sensor

**Technical Specification :**

Drop Counter sensor:

Use the Drop Counter to perform accurate, automatic titrations. Record the number of drops or calibrate the drop size to record volume of titrant added during a titration.

High accuracy, rubber sealed to protect from acids and bases

Measures up to a of 40 drops per second (even the very small drops less than 1mm)

LED drop indicator

Equipped with a UV filter that rejects all visible light



**Laboratory Instrument India**