



Email : sales@laboratoryinstrumentindia.com

Product Name :
Drop Counter Sensor

Product Code :
LBNY-0004-780055



Description :

Drop Counter Sensor

Technical Specification :

Drop Counter sensor:

Use the Drop Counter to perform accurate, automatic titrations. Record the number of drops or calibrate the drop size to record volume of titrant added during a titration.

High accuracy, rubber sealed to protect from acids and bases

Measures up to a of 40 drops per second (even the very small drops less than 1mm)

LED drop indicator

Equipped with a UV filter that rejects all visible light



Laboratory Instrument India