




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Redundant Truss Apparatus

**Product Code :**  
LBNY-0005-13900059



**Description :**

**Experiments:**

Study Of Strains, Stresses, Forces And Deflections In:

A Statically Indeterminate Structure

A Statically Determinate Structure

Operating Conditions

Operating Environment: Laboratory

Operating Relative Humidity Range:

80% At Temperatures Storage Temperature Range: 25c To +55c (When Packed For Transport)

Operating Temperature Range: +5c To +40c.

**Technical Specification :**

Redundant Truss Apparatus



**Laboratory Instrument India**