




Email : sales@laboratoryinstrumentindia.com

Product Name :
Redundant Truss Apparatus

Product Code :
LBNY-0005-13900059



Description :

Experiments:

Study Of Strains, Stresses, Forces And Deflections In:

A Statically Indeterminate Structure

A Statically Determinate Structure

Operating Conditions

Operating Environment: Laboratory

Operating Relative Humidity Range:

80% At Temperatures Storage Temperature Range: 25c To +55c (When Packed For Transport)

Operating Temperature Range: +5c To +40c.

Technical Specification :

Redundant Truss Apparatus



Laboratory Instrument India