




Email : sales@laboratoryinstrumentindia.com

Product Name :
Statics Work Panel

Product Code :
LBNY-0005-13900060



Description :

For use with Statics Fundamentals range, the Work Panel fits on most desk or bench tops. Fundamentals kits to the Work Panel to study or demonstrate one of the fundamental topics of static forces. Students or teachers fit the magnetic parts of their Statics.

Available Experiment Modules

Equilibrium Of Forces

Equilibrium Of A Beam

Suspension Cable Demonstration

Equilibrium Of A Rigid Body

Technical Specification :

Technical Operating Conditions

Operating temperature range: +5C to +40C

Operating relative humidity range: 80% at temperatures For use in: Well lit classroom or laboratory

Storage temperature range: 25C to +55C (when packed for transport).



Laboratory Instrument India

