



Email : sales@laboratoryinstrumentindia.com

Product Name :
Magnetic Stirrer with Heating

Product Code :
LBNY-0004-780056



Description :

Magnetic Stirrer with Heating

Technical Specification :

Magnetic stirrer with heating
Electronic speed control
White ceramic top plate
Housing and top plate resistant to acids, bases, solvents and scratches
Ergonomic control panel separated from the heating plate
Variable speed - 50 to around 1500 rpm
With rod
White top for clearly identifying solution color changes
Supplied with the magnetic turbulent
Heating up to 300 degrees Celcius.



Laboratory Instrument India

