



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Magnetic Stirrer with Heating

**Product Code :**  
LBNY-0004-780056



**Description :**

Magnetic Stirrer with Heating

**Technical Specification :**

Magnetic stirrer with heating  
Electronic speed control  
White ceramic top plate  
Housing and top plate resistant to acids, bases, solvents and scratches  
Ergonomic control panel separated from the heating plate  
Variable speed - 50 to around 1500 rpm  
With rod  
White top for clearly identifying solution color changes  
Supplied with the magnetic turbulent  
Heating up to 300 degrees Celcius.



**Laboratory Instrument India**

