



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vertical Stand Apparatus

Product Code :
LBNY-0005-14000014



Description :

The profile has extruded grooves along its length which enable the support of a wide variety of fixings and fastenings. A large, thick aluminium base plate supports a Vertical Stand Apparatus post at its centre, which itself is manufactured from extruded aluminium profile. These cross members are independently moved to suit the wall brackets of the apparatus. Three horizontal cross bars are supplied, which allow the secure fixing of wall mounted experiments from the range of products.

Technical Specification :

Levelling feet
Base plate 500mm x 500mm x 12mm nominal thickness
Cross member 420mm long
Vertical post 120mm wide x 1620mm long.



Laboratory Instrument India