



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Wireless Colorimeter Sensor and Turbidity

**Product Code :**  
LBNY-0004-780057



**Description :**

Wireless Colorimeter Sensor and Turbidity

**Technical Specification :**

Wireless Colorimeter Sensor and turbidity:  
Transmittance range: 0 -100%  
Wave lengths between 650nm, 600nm, 570nm, 550nm, 500nm and 450 nm  
Include cuvette rack with 12 slots  
Include ? 9 cuvettes  
Turbidity range: 0-400 NTU  
Accuracy: +5% NTU  
Include adapter for connection to PC or laptop.



**Laboratory Instrument India**