



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wireless Colorimeter Sensor and Turbidity

Product Code :
LBNY-0004-780057



Description :

Wireless Colorimeter Sensor and Turbidity

Technical Specification :

Wireless Colorimeter Sensor and turbidity:
Transmittance range: 0 -100%
Wave lengths between 650nm, 600nm, 570nm, 550nm, 500nm and 450 nm
Include cuvette rack with 12 slots
Include ? 9 cuvettes
Turbidity range: 0-400 NTU
Accuracy: +5% NTU
Include adapter for connection to PC or laptop.



Laboratory Instrument India