



Email : sales@laboratoryinstrumentindia.com

Product Name :
Voltage Sensor

Product Code :
LBNY-0004-780058



Description :

Voltage Sensor

Technical Specification :

Voltage Sensor :

- Measure the electrical energy generated by an electrochemical battery
- Measure the energy generated by a wind turbine or other renewable energy source
- Measure the voltage across components in a simple circuit
- Measure the induced voltage in a wire coil as a magnet passes through
- Includes at least 2 crocodile clip adapters
- Color coded wire cables (red & black)



Laboratory Instrument India