



Email : sales@laboratoryinstrumentindia.com

Product Name :
Whirling Of Shafts Apparatus

Product Code :
LBNY-0005-14000035



Description :

Experimental Capability:

Whirling of a shaft with concentrated masses
Whirling of a continuous shaft
Whirling of a long thin shaft supported at its ends.
Whirling of a shaft with intermediate nodes

Features:

Precision electronic tachometer
Three steel specimen shafts
Four variable position bearing housings provide nodal points
Infinitely variable control shaft speed
Three guard rings to limit whirling deflection.

Technical Specification :

Supplied complete with AC Control Unit, tachometer, set of specimen shafts and set of inertia discs and bushes. Comprising a bench mounted light alloy bed 1.8m in length and fitted with high speed DC motor.

Laboratory Instrument India

